


HEALQU™ Foam Dressing



HEALQU™ Foam Dressing is a three-dimensional polyurethane foam structure composed of a hydrophilic polymer which absorbs exudates and maintains a moist wound environment to promote healing. The foam material prevents the dressing from adhering to the wound or surrounding skin so it is painless to remove.

Benefits

- Especially suitable for heavily-exuding wounds
- EO sterilized
- A variety of sizes available
- Highly absorbent (more than 10 times its weight of fluid)
- Breathable
- Provides a moist environment for wound healing
- Comfortable and conformable

Indications

HEALQU™ Foam Dressing is indicated for heavily exuding wounds including:

- Venous ulcers (leg ulcers)
- II-IV Bedsores
- Diabetic foot ulcers
- Donor site wound
- Abrasion wound
- Operative incision

AVAILABLE SIZES:

Product	Size	Part number
HEALQU Non-Bordered Foam Dressing	2 x 2 x 0.2" (5 x 5 x 0.5cm)	9020701
HEALQU Non-Bordered Foam Dressing	6 x 6 x 0.2" (15 x 15 x 0.5cm)	9020702
HEALQU Non-Bordered Foam Dressing	4 x 5 x 0.2" (10 x 12.5 x 0.5cm)	9020703