



## Calcium Alginate Wound Dressing

### Advanced Wound Care

Our Calcium Alginate Dressing is made of natural fibers derived from seaweed which has been shown to have natural healing properties. The highly absorbent dressing quickly absorbs exudates and forms a bacterial resistant, gel-like covering over the wound, maintaining a moist environment for ideal wound healing. The gel-like surface also keeps the dressing from adhering to the wound to minimize pain on removal.



### Usage

This type of dressing is best used on moderately to heavily exuding partial to full thickness wounds including (but not limited to):

- Surface trauma and local hemostasis
- Post-surgical and donor site wounds
- Venous/arterial leg ulcer
- Cavity wounds
- Pressure sores
- Diabetic ulcers

Ref Number	Description	QTY
9020201	2x2" (0.5cm x 0.5cm)	10/100
9020203	4x5" (10cm x 12.5cm)	10/100
9020204	4x8" (10cm x 20cm)	5/50
9020205	1x12" (0.25cm x 30cm)	5/50



For more info: 1.844.4HEALQU | [www.healqu.com](http://www.healqu.com) | [info@healqu.com](mailto:info@healqu.com)