




# HEALQU™

The best in advanced wound care

## Bordered Silicone Foam Wound Dressing

### Advanced, Multi-layered Protection

HEALQU™ Bordered Silicone Foam Dressing consists of a gentle silicone wound contact layer, a super absorbent pad layer, and a vapor permeable and waterproof film. It is an ideal dressing for advanced wound exudate management. The multi-layer construction facilitates dynamic fluid management to provide an optimal moist wound environment, which promotes wound healing, minimizes trauma and reduces pain during removal of dressing.



### Usage

Our bordered Silicone Foam Dressing is used for a wide range of wounds such as:

- 1st and 2nd degree burns
- Leg and foot ulcers
- Traumatic wounds
- Surgical wounds
- Pressure ulcers
- Donor sites
- Skin tears



Ref #	Description	QTY
9020504	3x3" (7cm x 7cm)	5/50
9020501	4x4" (10cm x 10cm)	5/50
9020502	6x6" (15cm x 15cm)	5/50
9020503	7x7" (18cm x 18cm)	5/50

For more info: 1.844.4HEALQU | [www.healqu.com](http://www.healqu.com) | [info@healqu.com](mailto:info@healqu.com)