




HEALQU™

The best in advanced wound care

Hydrocolloid Wound Dressings

Breathable Protection

HEALQU™ Hydrocolloid Dressing absorbs exudates from the wound and is permeable to both oxygen and water vapor but impermeable to bacteria. The dressing forms a protective gel that provides a moist wound environment for ideal healing. The substrate gel also helps minimize tissue trauma and pain during dressing changes.



Usage

Our Hydrocolloid dressing are is used for a wide range of wounds such as:

- Skin graft donor site wounds
- 1st and 2nd degree burns
- Venous leg ulcers
- Arterial ulcers
- Pressure ulcers
- Abrasions

Ref Number	Description	QTY
9020801 (Standard)	2x2" (5cm x 5cm)	20/3200
9020802 (Standard)	4x4" (10cm x 10cm)	10/800



For more info: 1.844.4HEALQU | www.healqu.com | info@healqu.com

Manufactured for: HEALQU™ | Jersey City, NJ 07306 |
1.844.4HEALQU | www.healqu.com | info@healqu.com
Made in China
Doc #: 704-02009-01

HEALQU™