



**HEALQU™**  
The best in advanced wound care

# Super Absorbent Dressing

## The best in high-capacity exudate management

This dressing is a high-capacity, absorbent dressing particularly suited for heavily exuding wounds. The dressing consists of a breathable but waterproof, nonwoven backing, a super absorbent pad layer, and a hydro-penetrating nonwoven contact layer. This dressing locks exudates away and provides superior fluid management.



## Usage

For heavily exuding wounds including but not limited to:

- Partial and full-thickness wounds
- Post-operative wounds
- First and second degree burns
- Leg and foot ulcers
- Diabetic ulcers
- Pressure ulcers
- Surgical wounds

Ref Number	Description	QTY
9022007	5"x5" (12.5cm x 12.5cm)	10/10
9022011	6"x8" (15cm x 20cm)	10/10

For more info: 1.844.4HEALQU | [www.healqu.com](http://www.healqu.com) | [info@healqu.com](mailto:info@healqu.com)