



HEALQU®

The best in advanced wound care

Hydroconductive Wound Dressing

Layers of protection

Our Hydroconductive Wound Dressing is a sterile dressing for use on moderate to heavy exuding wounds. This technologically advanced dressing consists of a top layer, a super absorbent layer, and a conductive absorption layer. The three-layer construction not only facilitates fluid management, it also provides an optimal, moist wound environment to promote healing. Superior exudate absorption and non-adherent material helps to minimize the risk of maceration and damage to peri-wound skin.



Usage

For moderate to high exuding wounds including but not limited to:

- Non-infected diabetic foot ulcers
- Dehisced surgical wounds
- Pressure ulcers
- Leg ulcers
- Donor sites

Ref Number	HCPCS	Description	QTY
9021105	A6197	4"x5" (10cm x 12.5cm)	10/10



For more info: 1.844.4HEALQU | www.healqu.com | info@healqu.com